

SHOWER TROLLEY (Hydraulic)

Model	Size (LxDxH / mm)
HL-MT610	1900x640x(520~1000)
Angle	30 °
Weight	70kg
Load Capacity	120kg~150kg

- Patient and Assister feel comfortable when it is up and down by Hydraulic lift
- Regulate the angle of mattress make it easy to irrigation
- Easy to move patient with side mattress open
- Safety caster locking



SHOWER TROLLEY (Electric)

Model	Size (LxDxH / mm)
HL-MT620	1900x640x(520~1000)
Angle	30 °
Weight	70kg
Load Capacity	120kg~150kg
Power	AC 220V, 60HZ, 150W

- Bath equipment for the health of the patient
- Hydraulic Control of Height
- Easy to move patient with side mattress open
- Safety caster locking



물빠짐 배수로적용



물빠짐 배수관적용



환자 이동편리