

HL-SK-310

Hydraulic Patient Stretcher Cart



Side Rail Holder
Safety on patients get higher as the double locking system has been applied



Back Panel Angle Controller
The angle of back panel can be controlled within the range of 0° to 75° by the gas spring



Height Control Pedal
Height gets heightened when you continuously step on the "UP" pedal and gets down when you keep on stepping on the "DOWN" pedal for a while



Center Concentrated Brake
Each wheel has a center concentrated brake pedal, so that medical team can easily operate the brake at any position around the car



Training Wheels for Going Straight
Training wheels for going straight are installed in the middle of the car so that, it can be easily operated to go straight and also circumvolve on the spot treating the training wheels as its axis

Hydraulic Patient Stretcher Cart

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|--------|
| HL-SK-310 | 1940x666x500~835 | |

SPEC.

- Frame : Steel pipe
- Body : Aluminum
- Pump : Oil hydraulic pump
- Finish : Electrostatic powder coating(Use of the epoxy-polyester)
- Casters : 6inch swivel caster(Central locking system)
- Top : P.P
- Base : Steel pipe

Std. Acc.

- I.V. pole 1EA • Mattress 1EA

