

FOREST+LIFT

- Better social interaction
- Better access



- Adjustable to your needs
- Excellent outdoor and indoor performance
- Stability and easy to handle
- Quality & reliability
- Attractive design



- Better access



FOREST + LIFT

03

Excellent INDOOR and OUTDOOR performance for more independence and mobility.

- Excellent maneuverability
Motor 2 x 200W or 2 x 350W
- Excellent outdoor stability. The stability axle keeps the wheels in touch with the ground on uneven terrains up to 120 mm.
- Driving range: 30 km

01

300 mm lift for a better SOCIAL INTEGRATION

The new Forest Lift allows me to interact with people around me, enhancing my mobility and social integration.

02

300 mm lift for a better ACCESSIBILITY

- Kitchen cupboards
- Cash dispensers
- Bar
- etc...

04

ADJUSTABLE to your needs

- Armrests adjustable in height and width
- Adjustable seat depth
- Reclinable seat and back
- Adjustable controller
 - Leg rests adjustable in length and width
 - Foot rests adjustable in height

05

QUALITY & attractive design

- Robust quality of parts and frame supports efficient use in demanding outdoor and indoor situations.
- Attractive looks and quality design



Technical Information Forest lift

total width	seat width	seat/back angle	height armrests	height seat	total height	total length	seat depth	height back	users weight	weight
580	390 - 490	4° - 31°	200 - 250	530 - 830*	1130 - 1430	830**	420 - 480	510 / 610	130	155
650	440 - 540	4° - 31°	200 - 250	530 - 830*	1130 - 1430	830**	420 - 480	510 / 610	130	
650	500 - 600	4° - 31°	200 - 250	530 - 830*	1130 - 1430	830**	420 - 480	510 / 610	130	

Sizes in mm, weight in kg - * Without seat cushion - ** Without foot rests - Deviations in sizes and weight (15 mm, 1,5 kg, 1,5°) are possible.

batteries	motor	speed	driving range	inclination slope up	inclination slope down	obstacle height	turning radius	wheels
2x12V 70Ah	2 x 200 W or 2 x 350 W	6 km/h or 10 km/h	30 km	10°	10°	70	1400	front: air 260x85 rear: air 340x70

Class B wheelchair for in- and outdoor usage. - Not intended as a passenger seat during transport.



Opening up new horizons in mobility - www.vermeiren.com - info@vermeiren.com