

# Squod

## Stand UP



## Technical information Squod


										
total width	seat width	height armrests	seat height	total height	standing height	total length	seat depth	back height	max. users weight	weight
65	44	Min 20 & Max 28	58	107/117	143 148 163 168	40 45 50 55	50 50 60 60	130	134	

\* with backheight of 60 cm. Deviations in sizes and weight (1,5 cm and 1,5 kg) are possible



### Extra comfort

An ideal electronic wheelchair for inside and outside usage. The seat and back in Soft-version makes you feel comfortable and at ease in any situation. The collapsible controller comes in handy when you like to move closer to a table. This new Squod version offers a **complete stand-function** allowing you to raise yourself in all safety. Additional **knee support and a chest-safetybelt** will keep you steady at all times.



### The wheelchair must be adapted to you, not vice versa

Our basic principle is: the wheelchair must be adapted to you, not vice versa. This is why our wheelchairs have a modular design.

All types of wheelchairs can be adapted, have numerous accessories and all options are mutually interchangeable. This is the result of our continuous research and more than 51 years of experience.

## N.V. Vermeiren N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

Tel.: +32 (0) 3 620 20 20  
Fax: +32 (0) 3 666 48 94

[www.vermeiren.com](http://www.vermeiren.com) - [info@vermeiren.be](mailto:info@vermeiren.be)

For more information, contact your local dealer



**VERMEIREN**  
Opening up new horizons in mobility