











# Eclips XXL

SPECIAL BIG SIZES



## Technical information **Eclips XXL**

 total width	 folded width	 seat width	 height of armrests	 seat height	 total height	 total length	 seat depth	 back height	 weight
85	37,5	60	20	53	103	114	53	57	38,5
90	37,5	65	20	53	103	114	54	57	39
95	37,5	70	20	53	103	114	55	57	40
100	37,5	75	20	53	103	114	56	57	40,5

Deviations in sizes and weight (1,5 cm and 1,5 kg) are possible

### Special big sizes

Foldable wheelchair in aluminium with a double cross. The Eclips XXL brings you the original comfort with its anatomical seat and back. This wheelchair, with special reinforced rear wheels, can support a person up to 200 kg. The transfer of a person into the chair is very easy due to the removable armrests in aluminium.



*Brakes B74 and anti-tipping wheels B78 are optional*

### Your guarantee of quality

Regarding quality, the Eclips XXL complies with the highest norms laid down by certified test institutes in most European countries.

### The wheelchair must be adapted to you, not vice versa

Our basic principle is: the wheelchair must be adapted to you, not vice versa. This is why our wheelchairs have a modular design. All types of wheelchairs can be adapted, have numerous accessories and all options are mutually interchangeable. This is the result of our continuous research and more than 50 years of experience.

## N.V. Vermeiren N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

Tel.: +32 (0) 3 620 20 20  
Fax: +32 (0) 3 666 48 94

[www.vermeiren.com](http://www.vermeiren.com) - [info@vermeiren.be](mailto:info@vermeiren.be)

For more information, contact your local dealer



# VERMEIREN

Opening up new horizons in mobility

Subject to technical alterations