

FOREST + LIFT

Expand your limits



NEW

- Better social interaction
- Easy access to places difficult to reach

And also:

- Adjustable to your needs
- Excellent indoor and outdoor performance
- Stability and easy to handle
- Quality & reliability
- Attractive design
- Competitive price



FOREST + LIFT

Expand your limits

01

30 cm lift for a better SOCIAL INTEGRATION

The new Forest Lift allows me to interact with people around me, enhancing my mobility and social integration.

02

30 cm lift for a better ACCESSIBILITY

- Kitchen cupboards
- Cash dispensers
- Bar
- etc...



03

Excellent INDOOR and OUTDOOR performance for more independence and mobility.

- Excellent maneuverability
Motor 2 x 200W or 2 x 350W
- Excellent outdoor stability. The stability axle keeps the wheels in touch with the ground on uneven terrains up to 120 mm.
- Driving range: 30 km

04

AJUSTABLE to your needs












- Armrests adjustable in height and width
- Adjustable seat depth
- Reclinable seat and back
- Adjustable controller
- Leg rests adjustable in length and width
- Foot rests adjustable in height

05










QUALITY & attractive design

- Robust quality of parts and frame supports efficient use in demanding outdoor and indoor situations.
- Attractive looks and quality design

Technical Information Forest lift

 total width	 seat width	 seat/back angle	 height armrests	 height seat	 total height	 total length	 seat depth	 height back	 users weight	 weight
580	390 - 490	4° - 31°	200 - 250	530 - 830*	1130 - 1430	830**	420 - 480	510 / 610	130	155
650	440 - 540	4° - 31°	200 - 250	530 - 830*	1130 - 1430	830**	420 - 480	510 / 610	130	
650	500 - 600	4° - 31°	200 - 250	530 - 830*	1130 - 1430	830**	420 - 480	510 / 610	130	

Sizes in mm, weight in kg - * Without seat cushion - ** Without foot rests - Deviations in sizes and weight (15 mm, 1,5 kg, 1,5°) are possible.

 batteries	 motor	 speed	 driving range	 inclination slope up	 inclination slope down	 obstacle height	 turning radius	 wheels
2x12V 70Ah	2 x 200 W or 2 x 350 W	6 km/h or 10 km/h	30 km	10°	10°	70	1400	avant: air 260x85 arrières: 340x70

Class B wheelchair for in- and outdoor usage. - Not intended as a passenger seat during transport.

