



Sagitta

GO

every

DAY

be free

on it

Active aluminium wheelchair.
durability and reliability .



SAGITTA

EXTREME



Opening up new horizons in mobility

Dynamic and Activity. If you are the person whose life is full of any activity. If your requirements as for the wheelchairs are coming from the highest level-we have for you sagitta



...impossible becomes possible..

short



long



- **Balance** - perfect balance gives excellent manoeuvre.
- **Frame construction**- 'gold division' optimal proportions between the parts which are adjustable and solid gives extremely strong frame and beauty itself as well . Ideal order , like in cosmos, in the constellation SAGITTA.
- **Light and compact**- sagitta is made from aluminium -total weight 9,5 kg-it will be accompany you in every moment of yor life. that wheelchair was mady by support of the people which are the user of such a kind of the wheelchairs.
- **Camber** -you can choose between 0 , 2 ,4, 6, 8 degrees.
- **Back angle**-unique system of changing of the back angle in smooth way.
- **Seat and back**-tension changeable by velcro straps.
- **High of back**- changeable.
- **Upholstery** - made from CORDURA® : 10x more durable than cotton duck , 3x more durable than standard polyester ,2x more durable than standard nylon.

frame color:



| | | | | | | | | | |
|----------------------------------|--|----------------------------|---|----------------|--|---------------|----------------------------------|------------|-----------|
| | | | | | | | | | |
| 36 38 40 42 44 46 | 55-68 57-69 59-71 61-73 63-75 65-77 | 0° 2° 4° 6° 8° | 48 (without cushion) 52 (with cushion) | 55 60 65 | 32 34 37 39 42 44 46 | 10-poprzeczka | 23-28 28-33 33-38 38-43 | >9,5 kg | 100 kg |

Note: Specifications can change depending on the equipment and mode of wheelchair!



Opening up new horizons in mobility



back without push handles - standard



push handles - option



adjustable push handles - option



CORDURA®
Only by DuPont

standard



double footplate option



one piece aluminum plate SHORT - standard



aluminum open loop footrest - standard



aluminum cover plate - option



B28-standard cushion



Viscoelastic cushion - option



ROHO-option



back ERGO- Aluminium ergonomic back provides lateral stability-option



B78-anti-tips wheels option



Land cruiser tire Schwalbe-off road wheel 24x2-option



Right Run tire Schwalbe, wheels 22',24' - standard



wheels 24', 22' standard



wheel 24'x2 off road-option



wheel 24' exclusive -option



3'



4'



B85-transparent



B85-sagitta



B85- flowers 1



B85- flowers 2



B85- extreme



B85- pirat



standard push to lock -option



sport brake under the seat-standard



5'



Opening up new horizons in mobility