

ALK06-B03P-B

Three Function Electric hospital bed



Features:


- 1) Cold steel plate wholly molded surface, connected with arc soft join.
- 2) Streamline headboard and banister of PE.
- 3) 3pcs electro-motor, 1pc control box, 1pc handset
- 4) Automatic Long and Short lifting railing, safe and reliable, and can be fixed upward and downward
- 5) One set of handset
- 6) Four wheels adopt central control system, which is stable and reliable

Dimensions:

- 1) Size: L2260 × W1100 × H480-720mm
- 2) Angle of back section: 0 ~ 60° (± 5°)
- 3) Angle of leg section: 0 ~ 30° (± 5°)
- 4) Adjustable height: 480-720mm

Accessories:

- 1) I.V.pole
- 2) Wheels
- 3) Automatic lifting railing



3pcs electro-motor, 1pc control box, one set of handset, Battery back up

Can change the plug and voltage of the motor as inquiry



Two Sections I.V. Pole (Four Hooks)
Frame: Stainless steel

Twin-wheel Lockable Castors Central Brake System



Bed head and foot board
Frame: PE Molding
Movable head and foot boards with Lock-Switch, the board can be detached easily and quickly.



Accessories (option)

Mattress
Inside: high dense sponge, coco & palm fibre
Surface: waterproof & wearable cloth (fireproof as optional)



Automatic Long and Short lifting railing
Frame: PE molding, safe and reliable, and can be fixed upward and downward
Degree show function



Accessories (option)

Bedside table
Movable, can adjust the position low or high
Made of ABS and aluminum alloy