


# ALK06-A329P-B

## Luxury three crank hospital bed



**3**  
Cranks

### Features:


- 1) 3sets manual crank system. Cold steel plate wholly molded surface, connected with arc soft join.
- 2) Automatic lifting guardrail, safe and reliable, and can be fixed upward and downward
- 3) Four truckles adopt central control system, which is stable and reliable, 2pcs urine hooks, 6pcs IV drip holes
- 4) Integrated stamping frame and treated with electrostatic spraying, anti-aging and anti-rust
- 5) Max Loading: 280 Kg

### Dimensions:

- 1) Size: L2260 × W1100 × H480-720mm
- 2) Angle of back section: 0 ~ 60° ( ± 5° )
- 3) Angle of leg section: 0 ~ 30° ( ± 5° )
- 4) Adjustable height: 480-720mm

### Accessories:

- 1) I.V.pole
- 2) Wheels
- 3) Automatic lifting guardrail



Triple manual crank system  
The crank adopt Hidden design



Two Sections I.V. Pole (Four Hooks)  
Frame: Stainless steel  
Twin-wheel Lockable Castors Central Brake System



Bed head and foot board  
Frame: PE Molding  
Movable head and foot boards with Lock-Switch, the board can be detached easily and quickly.



**Accessories (option)**  
Mattress  
Inside: high dense sponge, coco & palm fibre  
Surface: waterproof & wearable cloth (fireproof as optional)



Automatic Long and Short lifting railing  
Frame: PE molding, safe and reliable, and can be fixed upward and downward  
Degree show function



**Accessories (option)**  
Bedside table  
Movable, can adjust the position low or high  
Made of ABS and aluminum alloy