

# ALK06-A235B

## Two crank hospital bed



### Features:



- 1) Double manual crank system
- 2) ABS bed head, 6 aluminum bars fence can go up and down.
- 3) 5" universal castors (2pcs with brake, 2pcs without brake)
- 4) 4pcs soft crashworthy bumpers, 2pcs urine hooks, 6pcs IV drip holes
- 5) The plate is made of stainless steel
- 6) Max Loading: 280 Kg

### Dimensions:

- 1) Size: L2130 × W980 × H500mm
- 2) Angle of back section: 0 ~ 70° (± 5°)
- 3) Angle of leg section: 0 ~ 30° (± 5°)

### Accessories:

- 1) I.V.pole
- 2) Wheels
- 3) Aluminum alloy protective railing
- 4) Oddment hook



# ALK06-A238P

## Two crank hospital bed

### Features:

ABS head and foot end, easy for installment.  
 Frame made of epoxy-coated Steel  
 Aluminum alloy protecting rails, easy and safe for movement up and down  
 The bed plate is made of grid.  
 Four castors, two with brake, two without brake  
 Weight capacity: 250 Kgs

### Dimensions:

- 1) Size: L2130 × W980 × H500mm
- 2) Angle of back section: 0 ~ 70° (± 5°)
- 3) Angle of leg section: 0 ~ 30° (± 5°)

### Accessories:

- 1) I.V.pole
- 2) Wheels
- 3) Aluminum alloy protective railing
- 4) Oddment hook

