

ALK06-B878C

Two crank manual Home care bed



Functions:


- 1) Double manual crank system
- 2) Solid wood outside, powder spray painting steel
- 3) Aluminum alloy bedside rails

Dimensions:

- 1) Size: L2080 X 1220 X 420mm
- 2) Angle of back section: 0 ~ 80° (± 5°)
- 3) Angle of leg section: 0 ~ 40° (± 5°)

Accessories:

- 1) I.V.pole
- 2) Aluminum alloy protective railing



Two Sections I.V. Pole (Four Hooks)
Frame: Stainless steel

5" universal castors(2pcs with brake,
2pcs without brake)



Collapsible side rail
6 aluminum bars fence with Lock-Switch
can be collapsed easily and quickly



Accessories (option)

Mattress

Inside: high dense sponge, coco & palm fibre
Surface: waterproof & wearable cloth
(fireproof as optional)



Accessories (option)

Bedside table

Movable, can adjust the position low or high, made of composite board and aluminum alloy, there is a rail on the composite board