



ALK06-AT1

Over bed table

Description:

- 1) ABS material injected plastic board
- 2) Aluminum alloy pull and stretch design

Size:

L900/1060 × 350mm



ALK06-AT2

Bedside table

Movable, can adjust the position low or high
Can be made as assembling
Specifications: 750 × 400 × 750-1050mm



ALK06-AT3

Bedside table

Movable, can adjust the position low or high
Made of ABS and aluminum alloy
Specifications: L800 × W400 × H750-1050mm



ALK06-AT5

Bedside table

Movable, can adjust the position low or high, made of composite board and aluminum alloy, there is a rail on the composite board
Size: 800*400*750/1050