

LS-EA5008A

Electric Bed



- 1** Specification
Overall size:2310mm(L) x 1190mm(W) x 510mm--700mm
Newly design of decent PP plastic head&foot board, rails.With easily operated membrane switch.
4 PCS 150mm Central-control system castors.
- 2** Function
Five functions.
This bed adopt luxurious membrane switch ,easily operated and pressed comfortably. With the "reset" button to make the bed be flat from any position soon . Also,a lock button with indicator light to show work status .
- 3** Hi--Low: 700mm-510mm
Back-rest lifting angle:0-65°
Knee-rest lifting angle:0-45°
Trendelenburg & Rev.Trendelenburg angle:0-10°
- 4** Loading Capacity: 250KGS
- 5** Standard Accessory: I.V. Pole
- Packing:2080mm x 1060mm x610mm

