

LS-EA5003B

Electric Bed



- 1 Specification**
Overall size:2170mm(L) x 950mm(W) x 480mm--680mm
Decent PP plastic head&foot board, rails.With easily operated membrane switch.
4 PCS 150mm Central-control system castors.
- 2 Function**
Five functions.
This bed adopt luxurious membrane switch ,easily operated and pressed comfortably. With the "reset" button to make the bed be flat from any position soon . Also,a lock button with indicator light to show work state .
- 3 Hi--Low : 680mm-480mm**
Back-rest lifting angle:0-65°
Knee-rest lifting angle:0-45°
Trendelenburg & Rev.Trendelenburg angle:0-10°
- 4 Loading Capacity: 250KGS**
- 5 Standard Accessory: I.V. Pole**
Packing:2080mm x 1060mm x 610mm

