

## LS-MA3050

Orthopedics traction Bed



- 1 LS920B-2 or LS920TW ABS Head / Foot board for choice
  - 2 LS-1500A or LS-1500B protective rails for choice
  - 3 Plastic and steel bed board mixed use
  - 4 6227-125mm diameter castor(2PCS), 6228-125mm diameter castor(2PCS)
  - 5 Stainless steel frame, easy to move and clean
  - 6 Left and right leg part position can be raised up by cranks, convenient for treatment and use
  - 7 Standars Accessories: I.V. Pole
- Weight: 125KGS

