



LS-3C

Manual Transfer Stretcher



- 1** It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;
- 2** Bright arrow symbols direct the users how to operate clearly.
- 3** Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easy control the stretcher direction
- 4** Hi-Low range: 860mm-570mm.
- 5** Standard accessories: mattress, I.V. Pole.
- 6** Packing: 0.98 CBM/PC.
- 7** Weight: 79KGS.

