

LS-4C

Hydraulic Transfer Stretcher



- 1** It is a type of hydraulic transfer stretcher. Five functions available by pressing the foot-pedals;
- 2** Bright arrow symbols direct the users how to operate clearly;
- 3** Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easy control the stretcher direction;
- 4** Hi-Low range: 825mm–580mm.
- 5** Standard accessories: mattress, I.V. Pole.
- Packing: 1.1CBM/PC.
- Weight: 82KGS.

