



LS-850E

Emergency Trolley

- 1 Material: Aluminum-Alloy, ABS plastic and stainless steel
 - 2 Overall size: 850mm x 520mm x 1010mm
 - 3 The whole trolley is including 5 layers of drawers, 2 small size (70mm in height), 2 medium size (140mm in height), and 1 large size (210mm in height). Drawer mainly consists of bottom panel board, frame and plastic medicine tray (size adjustable)
 - 4 Other components: 2 pcs of garbage bins, 1 pc of transparent file box, 1 pc of sliding side board, 1 pc of metal salver and 1 pc of I.V. Pole.
- Weight: 62KGS



LS-850Eb

Emergency Trolley



- 1 Material: Aluminium-Alloy, ABS plastic and stainless steel.
 - 2 Overall size: 850mm x 520 mm x 1010mm
 - 3 Newly matched with CPR board, Oxygen Cylinder holder, Power device, Needle disposal bottle.
 - 4 Components: 5 Layers of drawer, stainless top plate, 2 pcs garbage bins, 1 pc transparent file box, 1 pc sliding side board, 1 pc metal salver, 1 pc I.V. Pole.
- Packing: 910mm x 650mm x 1160mm
- Weight: 72KGS