



## LS-850M

Medicine Trolley

- 1 Material: Aluminum-Alloy, ABS plastic and stainless steel
  - 2 Overall size: 850mm x 520mm x 1010mm
  - 3 The whole trolley is including 6 layers of drawers, can be pulled out from both sides, self-locked when closed. Totally 50 medicine boxes on all drawers, on front of each medicine box, case-history card stucked.
  - 4 Other components: 2 pcs of garbage bins, 1 pc of transparent file box, 1 pc of sliding side board, and 1 set of five-hooks hanging device.
- Weight: 81KGS

Trolley