



LS-700D
Cure Trolley

- 1** 2 layers of big drug trays, you can put kinds of drugs or medical tools there; one big size drawer and two slide boards can be easily pull out and in;
 - 2** 1 I.V.Pole with 8 hooks, providing functional use of bottle and goods hanging;
 - 3** 4 white aluminum castors, 2 with brake, 2 without; the stainless-steel handrail makes the convenient movement of the trolley possible.
 - 4** Bright and clean color, makes the trolley looking very beautiful
- Oversize: 850mm x 650mm x 1100mm
 - weight: 33KGS



LS-800A
Medicine Trolley



- 1** The trolley consists of 6 layers, existing 5 cases of medicine boxes on both sides of each layer, totally 60 cases of the medicine boxes. At the bottom, there is a big waste-storage bin. Central-lock braking system, wheel diameter:125mm.
 - 2** Marble surface, so bright and easy to clean.
- Weight: 75KGS

