



LS-700A
42 Cases History Trolley



LS-700C
Medicine Trolley

1 The trolley consists of 42 pcs of case-history folder and 1 big drawer, all can be locked to help to keep the patient's history data inside safe. 4 pcs of 125mm castors, two with brake, and two without.

■ Weight: 65KGS.

1 The trolley consists of 6 layers of small drawers and 3 layers of big drawers, supplying various medicine boxes for different storage use. 4 pcs of 125mm castors, two with brake, and two without.

2 Standard accessories: 1 pc of garbage bin

■ Weight: 77KGS.

