

**EMERGENCY**

**TREATMENT UNIT**

*Quality and Service At It's Best*



**DF-900**

ISO 9000

ISO 13485

CE 0434



**DOCTOR'S FRIEND** medical instrument co., Ltd.

# EMERGENCY TREATMENT UNIT DF-900



**First-Aid Ambulance is the necessary appliance in modern versatile medical equipment. This instrument is equipped with multiple first-aid equipment, movable design and high mobility. It's the only choice at first-aid.**

## MULY-FUNCTIONAL INSTRUMENT TABLE :

- The top is equipped with the hanger which can hang mask, nose tube, etc. Oblong disinfection instrument box. Big round cotton pot. 2pcs of stainless steel instrument trays for storing emergency instruments. 4pcs of small cotton pots for loading gauze and cotton balls.
- Acrylic dust-proof cover which can prevent instruments from being polluted and is easy to take items.
- Halogen lamp DC-12V 35W in goose-neck tube design is adjustable at any angle with strong and concentrated light.

## MULTIFORM DEPOSITING VACANCY :

- 2 drawers with interior acrylic partition for conveniently depositing emergency medicines, syringes, needles, tablets, etc. Easy to handle, select and classify.
- 1 drawer for depositing the emergency instruments such as IV SET, funnel, respirator, suction catheter, gloves, etc. It's equipped with a lock to prevent expensive medicines & instruments from losing or be stolen.

## MOVABLE WRITING AND EMERGENCY BOARD :

- The movable & adjustable writing board can be used for filling in medical record or depositing heart beat monitor, electric shocker, etc.
- The movable heart emergency massage board can be used for C.P.R. Because the mattress of sick-bed and diagnosis bed of general hospital is too soft which can reduce emergency effect, the successful rate will increase if massage board is first put on the back and then proceed first-aid.
- The tong container installed at the side is easily removed to be disinfected.

## APPLIED RANGE:

Emergency room, recovery room, operating room, I.C.U., medical room, nurse station, caring center, etc.

## MATERIAL USED:

All made of specially treated stainless steel. Non-reflection, non-refraction, easy maintenance.

## STANDARD ACCESSORIES:

1. Multi-functional hanger ..... 1set
2. Oblong disinfection instrument box with cover ..... 1pc
3. Big round cotton pot ..... 2pc
4. Small round cotton pot (6cm) ..... 4pcs
5. Instrument tray (25cmx33cm) ..... 2pcs
6. Movable writing board ..... 1pc
7. Emergency board for CPR use .... 1pc
8. Big depositing box (at the back of machine) ..... 1pc
9. 500L oxygen cylinder ..... 2pcs
10. Pressure gauge, flow gauge, moisture vessel ..... 1set
11. Air way (S) ..... 1pc
12. Air way (L) ..... 1pc
13. Suction tube (180cm) ..... 1pc
14. Suction catheter ..... 1pc

## EMERGENCY OXYGEN UNIT :

- 2pcs of 500L oxygen high pressure steel cylinders is installed at the side of the emergency treatment unit. The 2 steel cylinders can be filled with oxygen by turn to efficiently bring the first aid into fully play. It's equipped with 1 set of pressure gauge, flow gauge, moisture vessel and oxygen mask.
- 3 sets of power sockets providing for monitor, electric shocker, etc. to prevent power not enough trouble while emergency.
- The depositing box installed at the back of the unit for storing suction catheter, X-ray photograph, etc.

## POWERFUL SUCTION UNIT :

- Equipped with super-powerful 200w oil-less vacuum pump, vacuum gauge and suction regulator. The volume of suction bottle reaches 2000c.c. It's installed with automatic full level control to secure life of motor.
- The depositing vacancy can store glucose, gauze, flexible bandage, cotton buds, etc.
- Anti-collision bumper bar can prevent collision while emergency transport and is equipped with wheel brake for fastening.

Power: AC 120 / 240V

Size: L 55 x W 45 x H98 (cm)

Weight: 55 KGS (ABOUT)



**DOCTOR'S FRIEND** medical instrument co., Ltd.