

# seca robusta 813

The standing,  
taking everything lightly



- Solid steel frame. For persons up to 200 kilograms.
- Non-slip rubber base. For safe standing during measuring.
- Big LCD figures. For easy reading.
- Step-off function. For weighing without switching on scales.

## seca robusta 813:

The standing, making everything look a breeze.



**T**ake a look at the extra-wide weighing platform. Even very big people can stand here safely and securely. People of larger girth too can be accommodated. Our intention is to lighten the task for weighty people.

A scale of standing, making the sturdy look beautiful.

More and more people weigh more and more. Now there's a stable scale for them. Thanks to its solid steel frame, it can withstand 200 kilograms and makes light of heavy use over many years. It's also no burden to step onto the scale, since its construction is extremely flat. That's a big advantage for people whose size makes it difficult for them to get about.

The standing, making the robust look elegant.

All rubber is not created equal. For this weighing platform seca has chosen a high-quality two-component nubby material which is soft like rubber and pleasant to the touch. But don't be fooled. This material is also tough. It fulfils two requirements at the same time: pleasant weighing and an easy-to-clean surface that stands up to all commercially available cleaning and disinfecting products.

The standing, paying attention to details.

Weight and girth are no problem on the seca robusta 813. You don't even have to bend over to switch it on, thanks to the unique step-off function. An additional problem with many scales is the display. Not so here, because thanks to the 29-millimetre-high LCD figures, it is extremely easy to read the measurements. You see, a scale from seca makes life easier by paying attention to details.

*The extra broad platform, 400 x 45 x 300 mm, and flat shape are weighty arguments, especially for heavy people.*



*High-quality, two-component nubby material feels as soft as rubber.*



*Big LCD figures. For easy reading.*

## seca robusta 813

### Technical Data

- Capacity: 200 kg / 440 lbs
- Division: 100 g / 0.2 lbs
- Dimensions (WxHxD):  
433 x 47 x 373 mm / 17 x 1.9 x 14.7"
- Weight: 2.6 kg / 6.3 lbs
- Functions: auto-HOLD, step-off, tip-on automatic switch-on, automatic switch-off
- Parts to be delivered: 4 AAA batteries, 4 extension feet for carpets