


UTILITY PLUS CART



- Designed to bring dressing set, instruments and medicine with easy when doctors and nurses for treating the patients
- Use stainless and ABS material for sterilization and proof infection
- Well stabilized and strong structured (use corner ring)
- Each drawer has ABS baskets with well organizing dividers
- Calmed and good rolling 4"casters with brake
- All stainless steel 304 construction with non spills edge and handle

SPECIFICATION

Body	1.2T SUS 304/27 plate
Shelf	1.2T SUS 304/27 plate
Drawer	ABS plastic & Aluminum
Basket	ABS plastic
Leg	Ø25.4x1.5T SUS 304/27 pipe
Guard rail	Ø9.6 SUS 304/27 bar
Cornering	Die casting
Bumper	Ø90 round rubber bumper
Caster	4 inch

Option Accessories



- Multi bin container
- IV holder
- Disposal Can
- Needle disposal holder



- Computer shelf
- Swivel laptop shelf.



- Disposal can
- Size : φ 120x140mm
- φ 200x210mm
- φ 250x260mm



- Oxygen tank & holder



- Shelf for waste disposal box