


ICU CART (Swing-out Type)



● SPECIFICATION

Top plate	Polypropylene
Body	Aluminum column & ABS plastic
Drawer	ABS plastic & Aluminum
Basket	ABS plastic
Handle	Powder coated steel
PC box	Powder coated steel
I.V pole	SUS
Side storage	ABS plastic
Dust box	ABS plastic
Top System I	Powder coated steel & ABS plastic & SUS
Caster	5 inch

- ICU Cart is designed for ICU (intensive care unit) which is Monitoring very il or very injured patients 24 hours and preparing first-aid care system anytime in emergency case.
- The small unit footprint maximizes storage, yet takes up less floor space. They side storage units conveniently swing out when needed, then close back into the cart reducing the footprint and providing a handle for maneuverability.

Model	Size (LxDxH / mm)	Remark
HL-ICU0908-1	755(1045)x520x1030	

● Standard Accessories



- Side storage (E type)
The view, only including tilt-bin in options



- Side storage (F type)
The view of including tilt-bin & 2 dust box in options



- Side storage (G type)
The view of including tilt-bin & 1 dust box in options



- Handle



- Multiple power outlet
220V Concent