

I.V Pole Stand

I.V Pole Stand with circle or rectangular or triangle handle enables patients to move conveniently and safely without more efforts for patient's treatment or private business. Constructed of stainless steel 304 pipe for rust proof and sanitary purpose. Durable materials and modern design to match hospital environment.

- All stainless steel 304 pipe $\Phi 25$ frame
- Telescoping tube & strong tubular stand, stainless steel
- I.V pole with four hooks, stainless steel
- Plastic cap to fix I.V hooks
- Two assistant hanger(SUS) coated rubber for 500cc, 1000cc bottles
- Screw knob for height adjustment, slow decent via inside buffer
- Height adjustable from 1300 to 2200 mm
- A.B.S plastic circle or triangle handle or A.B.S plastic rectangular handle with removable plastic tray
- Heavy-duty base of die casting five or four legs with 2" swivel casters
- Heavy-duty base of powder coated steel four legs with 2" swivel casters
- Available with non-magnetic I.V pole stand



I.V POLE STAND

Model	Size (LxDxH / mm)	Remark
HL-IPS6-0	$\Phi 700 \times (1300-2200)$	
HL-IPS6-N	$\Phi 700 \times (1300-2200)$	Non-magnetic

Option Accessories



● Needle disposal holder
HL-NDH2-AC



● IV pole holder (4Hooks)



● Tissue shelf
HL-TC245-AC



● Wire basket
HL-WB-AC



Disposal can
HL-DC122-AC



● Chest tube holder
HL-CTH-AC



● Utility holder
HL-UH96-AC



● Drain glass
HL-DG-AC



● Utility holder
HL-UH96-AC



● Oxygen tank holder
HL-OTH115-AC



I.V POLE STAND

Model	Size (LxDxH / mm)
HL-IPS4-0	$\Phi 500 \times (1300-2200)$



I.V POLE STAND

Model	Size (LxDxH / mm)
HL-IPS4-3	$\Phi 600 \times (1300-2200)$

- Not movable