



ENDOSCOPIC CABINET (Holder Type)

Model	Size (LxDxH / mm)	Remark
HL-ECK-M	900x450x2100	6 Holder
HL-ECK-L	1200x450x2100	8 Holder



▲ heater & fan



▲ Endoscope holder (Telescopic runner type)



▲ Endoscope holder (Panel type)



▲ Catheter tube holder



ENDOSCOPIC CABINET (Panel Type)

Model	Size (LxDxH / mm)	Remark
HL-ECK-4	900x450x2100	4 Holder
HL-ECK-5	1200x450x2100	6 Holder

- Endoscopic Cabinet is designed for drying and sterilizing the cleaned endoscope by fan heater and U.V lamp.
- This is designed to keep high-priced endoscope safely and efficiently in examination center and endoscope examination room and provide more convenience with user.
- Heater is equipped in interior to dry moisture after cleaning endoscope.
- U.V lamp is equipped to sterilize to left and right sides and Timer to keep uniform temperature in chamber.
- All function operate automatically via sensor attached in sliding door when sliding door open and close.
- Locking unit in doors for user's convenience and adjustable foot in floor
- Heavy duty stainless steel 304 construction
- Stainless steel leg and adjustable urethane foot
- 5T strong clear glass and stainless steel door frame
- U.V lamp 20W (life time 5,000 hr)
- Power: AC220V 50/60Hz.
- Electric code: 2m
- Control box: On&Off power switch, Timer(30min.)
Heater & Fan(1Kw/hr), Counter
Air circulation fan, Main switch
Limit sensor, Hour meter(U.V lamp working time check)