


EMERGENCY CART



- Designed to manage a first-aid patient in emergency
 - Ergonomic designed push-handle and efficient structure
 - Total locking key system
 - Colors: Red, Yellow, Blue, Green
- From A-type to D-type, side storage has various kind of model, so you can choice a best adaptable type for your item
- Handy tilt-bins manufactured for fitting to side storage items can store various kind of nurse's items, more safely and usefully keeping all necessary items, then also easily exchanged when empty

● SPECIFICATION

Top plate	Polypropylene & Stainless steel
Body	Aluminum column & ABS plastic
Drawer	ABS plastic & Aluminum
Basket	ABS plastic
Sliding side table	ABS plastic
Guard rail	Aluminum
Caster	5 inch
Side storage	ABS plastic
Dust box	ABS plastic

● Standard Accessories



- Side storage (A type)
The view, only including tilt-bin in options



- Side storage (B type)
The view of including tilt-bin & dust box in options



- Push handle



- C.P.R board
(ABS Plastic)



- Sliding side shelf
The sliding side shelves shall be used to support small objects and document