

# Sunny & 925

Wheelchairs with a heart for children



Vermeiren has developed a separate series of wheelchairs specially designed for children. It is clear that these are more than simply smaller versions of adult wheelchairs. It is our aim to give the child as much freedom of movement as possible and to optimise the contact between the child, parents and therapists.

## Sunny

The electronic wheelchair for children which gives even greater freedom of movement

### Electronic and made-to-measure

The sunny is an electronic wheelchair for children which has a **microprocessor**. It is easy to operate, with a minimum of effort and maximum safety. Due to the **5 programs**, the wheelchair adapts itself to the specific needs of the child.

## 925

The trusty standard












### the companion for your child

The 925 is the number one in wheelchairs for children. It is the companion for your child and provides you with the assurance of quality. For example, the robust metal tubing can resist small knocks and special

clothing protectors ensure that clothing does not become entangled in the wheels of the wheelchair.



## Technical information Sunny & 925

										
total width	folded width	seat width	height of armrests	seat height	total height	total length	seat depth	back height	distance footrests-seat	weight
<b>Sunny</b>										
49	40	32	22,5/15	53	98	80	30/37	49	21,5/29	61 (batt.incl.)
51	40	35	22,5/15	53	98	80	30/37	49	21,5/29	61
54	40	37	22,5/15	53	98	80	30/37	49	21,5/29	61
<b>925</b>										
48	32	29	13/15,5	49	91	77	33	35	adjustable from 10 to 30 cm	13
54	32	35	13/15,5	49	91	77	35	35		13

Deviations in sizes and weight ( $\pm 1,5$  cm and  $\pm 1,5$  kg) are possible

## Sunny





### Your ideal travelling partner

The batteries, the back and the seat cushion can be **removed very easily**. As a result, the Sunny can be folded, it becomes considerably lighter and it can be easily stored away in the boot of the car.

### Sunny grows with your child

Children grow unbelievably fast. The Sunny adapts itself to your child: the depth of the **back**, the height of the **foot-piece**, and both the height and depth of the **armrests** can all be adjusted.

### Choice of 2 colour combinations

Frame	Upholstery
 yellow	 black
 red	 black

## 925

### Foot-plate can be adjusted up to 10°

Some children grow quicker than others. It is therefore important that the **foot-plate** can be adjusted in both **height** and in **depth**.




### Facile à transporter

The 925 is easy to take with you. It is light (only 13 kg) and can be folded. It is hence easy to take this wheelchair with you **in the car or the train**.

### Options

- adjustable legrests
- anatomical seat and back
- Quick release wheels

### Choice of 2 colour combinations

Frame	Upholstery
 yellow	 black
 light blue	 blue

## N.V. Vermeiren N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout, Belgium

Tel.: +32 (0)3 620 20 20

Fax: +32 (0)3 666 48 94

Visit our website: [www.vermeiren.com](http://www.vermeiren.com)



# VERMEIREN

Opening up new horizons in mobility

For more information, contact your local dealer.

