

HL-SK-153

Pediatric Hospital Bed ✓



Dining table

P.P materials integral table




Side rail

Prevention of a hurting from a falling by four(4) side rails of narrow vertical bar



Motor

Durability-tested 50,000 times LINAK actuator systems (made in Denmark)



Upper Panel Cover

Hinge and pin that connect the upper panel are covered so that the lifespan of mattress can be lengthened and your hands can be prevented from getting hurt while you replacing the mattress cover.



Pediatric Hospital Bed

Model	Size (LxDxH / mm)	Remark
HL-SK-153	2176x1020x385~715	

SPEC.

- Bed frame : Steel plate
- Bed top : P.P
- Head, Foot panel : P.P
- Finish : Electrostatic powder coating
(Use of the epoxy-polyester)
- Casters : 5inch swivel caster
- Motor : LINAK Actuator 3EA,
Control box 1EA, Hand set 1EA

Std. Acc.

- Over bed table 1EA
- Side rail 4EA
- I.V. pole 1EA
- Mattress 1EA

