




Side rail

Clamp fixed to a horizontal frame Excellent strength in both directions Sleek exterior side rails



Head/Foot panel

One-Touch Detachable Head/Foot panel




Dining table

P.P materials integral table




Tops

Hinge and Pin wrapped in a cover Extend the life of Mattress No risk of injury when changing Mattress cover





Motor

Durability-tested 50,000 times LINAK actuator systems (made in Denmark)



Nurse Control Panel (Option)



Electric Hospital Bed

Model	Size (LxDxH / mm)	Remark
HL-SK-151	2176x926x385~715	

SPEC.

- Bed frame : Steel plate
- Bed top : P.P
- Head, Foot panel : P.P
- Finish : Electrostatic powder coating (Use of the epoxy-polyester)
- Casters : 5inch swivel caster
- Motor : LINAK Actuator 3EA, Control box 1EA, Hand set 1EA

Std. Acc.

- Over bed table 1EA
- Side rail 2EA
- I.V. pole 1EA
- Mattress 1EA

