

# Streeter

*Progressive movement*



## Technical information **Streeter**

Red (S)	46-57	74	20	53	34	4,7	50
Beige (M)	47-60	90	17	55	36	5	65
Green (L)	60-73	96	21	59	40	5,7	80
Blue (XL)	72-90	102	21	63	43	6,4	90

Deviations in sizes and weight (1,5 cm and 1,5 kg) are possible.

### **Streeter** : *Flexible mobility*

The Vermeiren Streeter was designed to ease the learning and relearning of walking. The locking system on the backwheels will prevent the Streeter from moving backwards and will encourage the user to move forward.

The system is fully reversible so that your Streeter can be used in both directions. The aluminium frame is fully collapsible which makes it excellent for transport and storage. 4 different colors represent the 4 different sizes.

### **Your guarantee of quality**

Regarding quality, the Streeter complies with the highest norms laid down by certified test institutes in most European countries.

### **How to pick your Streeter**

The different models :

Red (S) - Beige (M) - Green (L) - Blue (XL)

It suffices to pick a Streeter where the height of the handles is equal to the height of the hips of the user. Every Streeter is adjustable in height for a perfect personal solution.



S



M



L



XL

## **N.V. Vermeiren N.V.**

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

Tel.: +32 (0) 3 620 20 20  
Fax: +32 (0) 3 666 48 94

[www.vermeiren.com](http://www.vermeiren.com) - [info@vermeiren.be](mailto:info@vermeiren.be)

For more information, contact your local dealer



# **VERMEIREN**

Opening up new horizons in mobility

Subject to technical alterations