











VERMEIREN

708
delight
vermeiren

A new standard



Technical information 708 Delight

 total width	 folded width	 seat width	 height of armrests	 seat height	 total height	 total length	 seat depth	 back height	 weight
57	31	39	22	51	92	110	43	43	16,5
60	31	42	22	51	92	110	43	43	16,5
62	31	44	22	51	92	110	43	43	17,5
64	31	46	22	51	92	110	43	43	17,5
66	31	48	22	51	92	110	43	43	17,5
68	31	50	22	51	92	110	43	43	17,5

Deviations in sizes and weight ($\pm 1,5$ cm and $\pm 1,5$ kg) are possible

The "new" standard

The 708 delight has removable or swing-back armrest. Besides that, the basic version is also equipped with frontwheel fork adjustment. The front wheels can be adjusted in two different positions. The steel was replaced by aluminium which results in a significant loss of weight.











Basic Features

- Removable armrests (Bo2 or Bo3)
- Removable legrests (Bo6) with calf strap (B57)
- Frontwheel fork adjustable in 2 positions
- Metal- or Plastic wheels with spokes (Rear wheels)
- Rear wheels 24" (Air or PU)
- Front wheels 8" (Air oder PU)
- Flexible seat and back

Options

- Armsrests adjustable in height (Bo5)
- Legrests adjustable in height (Bo8 or BZ7)
- Pushing model 12"-Rear wheels (T30)
- Seat width 50 cm
- Lowered seatingheight (SI)
- Firm, padded seat/back (L14/15 • L24/25)
- Brakes (B74)
- Brakes + removable wheels (B74 B80)
- Extensive options (Modular system)

Frame	Upholstery
 Blue C27	 Black
 Red C32	 Black
 Grey C29	 Black
 Green C26	 Black



For more information, contact your local dealer

N.V. Vermeiren N.V.

Vermeirenplein 1-15

B-2920 Kalmthout

Belgium

Tel.: +31 (0) 3 620 20 20 • Fax: +31 (0) 3666 48 94

Also visit our web-site: www.vermeiren.com



VERMEIREN

Opening up new horizons in mobility