

ALK06-A328L

Extra low three crank hospital bed



Features:


- 1) 3sets manual crank system
- 2) ABS bed head, 6 aluminum bars fence can go up and down.
- 3) 4pcs soft crashworthy bumpers
- 4) accessories: 5" universal castors (2pcs with brake, 2pcs without brake), 2pcs urine hooks, 6pcs IV drip holes
- 5) Integrated stamping frame and treated with electrostatic spraying, anti-aging and anti-rust
- 6) Max Loading: 280 Kgs

Dimensions:

- 1) Size: L2130 × W980 × H320-620mm
- 2) Angle of back section: 0 ~ 70° (± 5°)
- 3) Angle of leg section: 0 ~ 30° (± 5°)
- 4) Adjustable height: 320-620mm

Accessories:

- 1) I.V. pole
- 2) Wheels
- 3) Aluminum alloy protective railing



Triple manual crank system

The crank adopt Hidden design



Two Sections I.V. Pole (Four Hooks)
Frame: Stainless steel

5" universal castors (2pcs with brake, 2pcs without brake)



Bed head and foot board
Frame: ABS Molding
Movable head and foot boards with Lock-Switch, the board can be detached easily and quickly.




Accessories (option)

Mattress
Inside: high dense sponge, coco & palm fibre
Surface: waterproof & wearable cloth (fireproof as optional)



Collapsible side rail
6 aluminum bars fence with Lock-Switch can be collapsed easily and quickly



Accessories (option)

Telescopic Dining Table
The feature can shorten the length of table while not in use to provide against possible harms of protrusion