

# LS-1C

Connecting Transfer  
Stretcher For Operation Room

- 1 It is a manual connecting Transfer Stretcher, can be used to defend cross-infection;
- 2 Cart A inside the operation room, Cart B outside, adjusting them to the same height, then push Cart B to Cart A, two carts will be connected, following the orbit, bedboard C can be pushed from one cart to another, and VICE VERSA;
- 3 Bright arrow symbols direct the users how to operate clearly;
- 4 Hi-Low range: 920mm-645mm;
- 5 Standard accessories: mattress, I.V.pole;
- Weight: 145KGS

